

TS-700

TRANSPORT AND EMERGENCY STRETCHER

Ref.44EU2609



mubi
M E D I C A L

TS-700

TRANSPORT AND EMERGENCY STRETCHER



COMFORT AND DYNAMISM

TS-700 protects your patients while protecting your investment, thanks to a transport, examination and recovery stretcher that meets regulatory standards, ensuring quality care.

DESIGN AND SAFETY

TS-700

TRANSPORT AND EMERGENCY STRETCHER

Sliding backrest

reduces patient migration resulting in less boosting, eases burden on caregiver's back and helps reduce injuries, with a four-section platform: one fixed and three movable sections, removable and made of high-strength polymer.



Fold-away push handles

facilitate transport, improve maneuverability and provide access to the patient.



Knee gatch

provides the patient added comfort and prevents patient migration, minimizing boosting.



Bumper wheels

allows protect the structure and building walls institutions from possible impacts.

TS-700

TRANSPORT AND EMERGENCY STRETCHER



COMFORT AND SAFETY

Ideal for hospital environments, high strength steel frame, foam and laminate surfaces; a precise combination of materials for support, comfort and reliability.



Low height position
facilitates a safer patient
ingress and egress
reducing falls risk.



Recovery chair position
improves patient comfort, and
facilitates treatments (airway clearance
or administering medication).

TS-700

TRANSPORT AND EMERGENCY STRETCHER

DESIGN AND SAFETY

High resistance materials and easy disinfection finishes, compatible with products used in hospitals such as: hydrogen peroxide, quaternary ammonium, sodium hypochlorite and antiseptic alcohol**.

SAFETY SIDE RAILS



One-handed siderail release simple and efficient locking mechanism which avoids pinch point risks.

One-section, foldable, made of high-strength polymer and aluminum cover, easy to clean and disinfect, simple one-handed operation, facilitating medical care.

RETRACTABLE FIFTH WHEEL

allows optimal handling and control during transport.



CENTRAL-LOCKING BRAKE SYSTEM

prevent wheels from rolling or swiveling as safe patient ingress-exit.



** Recommended composition:

- Hydrogen peroxide 3% and 7.5% 200 ml peroxide per liter of water.
- Quaternary ammonium 4 ml per liter of water.
- Sodium hypochlorite 4 spoonfuls for each liter of water approximately.
- Antiseptic alcohol 60-80% 7 parts of alcohol x 3 parts of water.

TS-700

TRANSPORT AND EMERGENCY STRETCHER



Permanently mounted IV pole
at head-end, 2 hooks.



Removable IV pole
2 stage and 2 hooks.

TELESCOPIC IV POLE

One fixed and folding (left of backrest) and one removable to be placed in any of the 4 corners.

TS-700

TRANSPORT AND EMERGENCY STRETCHER

IDEAL FOR EMERGENCY
AND TRAUMA UNITS



Backrest



Horizontal



Trendelenburg



Reverse trend.



Knee gatch



Recovery Chair

*Ver cuadro de especificaciones.

STANDARD DESIGN WITH MULTIPLE CONFIGURATION
OPTIONS, WITH MOBILITY AND SAFETY IN MIND

TS-700

TRANSPORT AND EMERGENCY STRETCHER



Head of bed angle indicator
aids with airway clearance head elevation compliance.



Hydraulic elevation system
exceptionally durable and easy-to-clean anodised aluminium columns protect the mechanism.



Base cover
for use as an integrated oxygen tank holder or as a space for patient belongings.

TS-700

TRANSPORT AND EMERGENCY STRETCHER

SURFACE PAD (Optional with additional charge)

Helps in the treatment and prevention of pressure injuries, high density polyurethane (PU) foam core, $40 \pm 2 \text{ Kg/m}^3$ (3107.3 lb/ft³)*.



PAD COVER COLOR

TURQUOISE
HOSPITAL EXTREME

HOSPITAL EXTREME COVER

Made of nylon, permablok lacquer coating, light protection AATCC TM 16-3 / NTC 1479 / CFFA-2 a1 - Blue Wool, waterproof, anti-static ASTM D-257, bacteria resistance AATCC 147, JIS Z 2801, mildew resistant NO GROWTH (NG), invisible waterproofed seams, stitch resistant DIN 54301, flame retardant FMVSS 302 and latex free.

TS-700

TRANSPORT AND EMERGENCY STRETCHER



OVERALL DIMENSIONS *

Length	87.3 in (2210 mm)
Width	36 in (915 mm)

PATIENT SURFACE *

Length	77 in (1950 mm)
Width	28 in (710 mm)

MOTIONS *

Backrest	0° to 90°
Footrest	0° to 45°
Trendelenburg	17°
Reverse trendelenburg	-17°
Minimum height	23.2 in (590 mm)
Maximum height	34 in (863 mm)

* Motions performed by pneumatic and hydraulic systems, assisted by pedals on both sides.

** Tolerances: degrees $\pm 3^\circ$ dimensions 0.5 in (12.7 mm).

PACKAGING

- Container: 40HC – 39 units, 40STD – 26 units, 20STD – 12 units.
- Packing: Carton box 14795, Stowage 14796.
- Packing dimensions: Length 2290 mm Width 1020 mm Height 870 mm.
- Product weight + packaging 439.82 LBS (199.5 kg).

TS-700

TRANSPORT AND EMERGENCY STRETCHER

STANDARD FEATURES *

- 700 lbs weight capacity
- 28-inch wide patient deck
- 4-section polymer patient surface
- Retractable fifth wheel
- Knee gatch / recovery position
- Fold-away push handles on head-end
- Sliding backrest
- Permanently mounted ivpole at head-end, 2 hooks
- Angle indicators for head of bed and trend
- 8-inch integral casters with wheel covers
- Pneumatic assisted backrest
- Trendelenburg and reverse trendelenburg
- Central-locking brake system
- Brake and steer pedals on head and foot ends
- Hydraulic elevation control pedals on both sides
- Fold-down siderails
- 4 integrated IV pole sockets
- Drainage bag holders
- Bumper wheels on each corner
- Patient storage tray on base
- Oxygen tank holder on base

MANUFACTURING AND FINISHING

- Structure made of high-strength cold rolled steel tubing, reinforced with structural tubing.
- Finished with electrostatically applied, rust-resistant epoxy-polyester powder paint.

OPTIONAL ACCESSORIES (ADDITIONAL COST)

- 3-inch surface pad with antibacterial and anti-static cover
- 4-inch surface pad with antibacterial and anti-static cover
- Removable IV pole, 2 stage and 2 hooks

IMPORTANT NOTES

- ACCESSORIES Please consult with your executive.
- MUBI Medical reserves the right to change specifications without prior notice.
- Three-year Limited Product Warranty.

TRANSFORMING SPACES, IMPROVING LIVES



Medellín – Colombia
Calle 2 # 50-185
(57) 604 285 2838

Cali – Colombia
Calle 18 # 118-85
(57) 602 489 6999

